

TOC 1500

Air sans CO2

Numéro de pièce: 42-1030




La description

Conçu pour fournir une alimentation en air fiable aux analyseurs de COT, le TOC1500HP est un système robuste utilisé dans les laboratoires pour produire du CO2 propre et sec et sans CH4



Types d'application

Analyseurs de COT



Generator Details

- **Types de gaz:** Air
- **Compresseur intégré:** Non

Principales caractéristiques

- Produit un air propre, sec, sans CO2 et sans CH4
- Idéal pour les analyseurs de COT NPOC
- Fonctionne 24 heures sur 24, 7 jours sur 7, si nécessaire.
Montage mural pour les petits laboratoires
- Peu de pièces mobiles signifie peu de maintenance et garantit une longue durée de vie du générateur
- Fonctionnement totalement silencieux
- Configuration minimale requise
- Garantie complète sur site de 12 mois

Spécifications techniques

TOC 1500 Zero Air sans Co2

| | |
|-------------------------------|------------------------|
| Type de gaz | Air |
| Débit de gaz max: | 1500cc/min |
| Pression de sortie max: | 80psi/5,5bar |
| Pression de sortie max: | 80psi/5,5bar |
| Pureté max: | <1 ppm CO2/ 2.0ppm THC |
| Pression d'entrée min: | 90psi/6,2bar |
| Pression d'entrée max: | 145psi/10bar |
| Exigences de pureté de l'eau: | N/A |
| Consommation d'eau: | TBC |
| Voltage: | 100-230 ± 10% VAC |
| Fréquence: | 50 / 60 Hz |
| Courant: | 8 Amps |
| Accréditations: | CE |
| Taille (HxLxP) mm | 660 x 250 x 160 mm |
| Taille (HxLxP) pouces | 26 x 9,8 x 6,3 inches |
| Poids du générateur | 12kg / 26,4lbs |

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

