

# Solution totale Genius 8080

Numéro de pièce: XX-XXXX



## La description

Les Genius 1053 et Genius 3070 sont construits autour de la technologie éprouvée Peaks, fournissant une solution de gaz totale pour le triple quadripôle Shimadzu LCMS-8080.



## Generator Details

- **Types de gaz:** Azote
- **Compresseur intégré:** Non

## Principales caractéristiques

- **Flexible** - Disponible avec ou sans compresseur d'air interne
- **Pureté** - Zéro air avec moins de 0,1 ppm d'hydrocarbures (sous forme de méthane)
- **Durable** - Peu de pièces mobiles avec un minimum de maintenance
- **Solution à source unique** - Un générateur peut répondre aux exigences de plusieurs applications
- **Installation simple** - Générateur conçu comme plug & amp; système de jeu
- **Économique** - Plus rentable que toute autre méthode d'approvisionnement en gaz
- **Pratique** - Gaz à la demande, pas de danger pour la santé, pas besoin de s'inquiéter de manquer d'essence

## Spécifications techniques

### Solution totale Genius 8080

Type de gaz	Azote
Temps de mise en route:	60 mins
Courant:	0
Niveau de bruit:	0

# [PEAK Protected]<sup>TM</sup>

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: [www.peakscientific.com/ordering/](http://www.peakscientific.com/ordering/)
- For service plans visit: [www.peakscientific.com/service/service-plans/](http://www.peakscientific.com/service/service-plans/)

**Peak Scientific UK**  
Tel: +44 (0)141 812 8100  
Fax: +44 (0)141 812 8200

**Peak Scientific North America**  
Tel: +1-800-767-6532  
Fax: +1-978-608-9503

**Peak Scientific China**  
Tel: +86 21 5079 1190  
Fax: +86 21 5079 1191  
For a full list of our worldwide office locations, please visit:  
Web: [www.peakscientific.com](http://www.peakscientific.com)  
Email: [discover@peakscientific.com](mailto:discover@peakscientific.com)

Peak Scientific's Quality Management System conforms to: ISO:9001

